**Monoclonal Mouse
Anti-Human CD68**
Clone PG-M1
Code M0876

ENGLISH

Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD68, Clone PG-M1, is intended for use in immunohistochemistry (IHC). The antibody labels macrophages and is a useful aid for the classification of acute myeloid leukemia (AML), and histiocytic sarcoma (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.
Summary and explanation	CD68 is a highly glycosylated lysosomal membrane protein with an Mr of 110 000. The CD68 protein belongs to a family of lysosomal glycoprotein (LGP)/plasma membrane shuttling proteins that play a role in endocytosis and/or lysosomal trafficking. CD68 is expressed strongly in cytoplasmic granules, and weakly on the surface of macrophages, monocytes, neutrophils, basophils and NK-cells. Additionally, CD68 is expressed by some peripheral blood B cells and may be weakly expressed in B-cell type acute lymphoblastic leukemia (B-ALL) cells. CD68 can also be found in the cytoplasm of non-hematopoietic tissues, especially the liver, and renal glomeruli and tubules (2). Unlike many other CD leucocyte antigens, the CD68 molecule is antigenically very heterogenous, and different antibodies to CD68 show different cellular reactivities (3). Refer to <i>Dako</i> General Instructions for Immunohistochemical Staining or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> PG-M1 (3). <u>Isotype:</u> IgG3, kappa. <u>Mouse IgG concentration:</u> see label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
Immunogen	Human mononuclear spleen cell preparation containing more than 80% Gaucher's cells (3).
Specificity	The antibody was clustered as anti-CD68 at the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Boston in 1993 (4). The antibody labels COS-1 and WOP cells transfected with CD68 cDNA. Unlike other CD68 antibodies, which label both macrophages and myeloid cells, the PG-M1 antibody detects a fixative-resistant epitope on the macrophage-restricted form of the CD68 antigen (3).
Precautions	1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, Bouin's fixative or B5 fixative (3). Pre-treatment of deparaffinized tissues with chymotrypsin, trypsin, pepsin, pronase (3), proteinase K, or heat-induced epitope retrieval is required. For heat-induced epitope retrieval of tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone fixed, frozen sections, although the labeling is slightly weaker than that observed in paraffin sections (3, 4), and cell preparations (3). The user must validate the staining procedure.
Staining procedure	These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms. <u>Dilution:</u> Monoclonal Mouse Anti-Human CD68, Code M0876, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. <u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. <u>Visualization:</u> Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Product-specific limitations	The PG-M1 antibody labels some non-hematopoietic malignancies, especially about 10% of melanomas (3), therefore the classification of histiocytic sarcoma should always be supported by labeling of neoplastic cells for at least another macrophage-restricted marker and/or the leucocyte common antigen, and non-labeling of epithelial and melanoma-associated antigens (1).
Staining interpretation	Cells of the monocyte/macrophage lineage labeled by the antibody show a cytoplasmic (diffuse or granular) staining pattern.
Performance characteristics	<u>Normal tissues:</u> In normal peripheral blood only monocytes are labeled by the antibody. In a wide range of tissues tested, all macrophages were labeled. In bone marrow also osteoclasts were labeled, whereas granulocytes and myeloid precursors were consistently not labeled. A weak reactivity of some megakaryocytes was observed in about 20% of the cases. Kupffer's cells in the liver, mast cells, and synovial cells were the only additional normal cells labeled by the antibody (3). <u>Abnormal tissues:</u> Among 431 malignancies of the lymphohematopoietic system, reactivity with the antibody was restricted to acute myeloid leukemias of the M4 and M5 type, true histiocytic sarcomas, and mastocytosis. Consistently unreactive were acute myeloblastic leukemias of M1, M2, and M3 type, malignant non-Hodgkin's lymphomas of B- and T-cell type, Hodgkins lymphoma, acute lymphoblastic leukemias, and chronic myeloid leukemia. The majority of 370 non-hematopoietic tumors were not labeled with the antibody, exceptions being 15/15 granular cell myoblastomas, 6/13 kidney clear cell carcinomas, 4/10 glioblastomas, 10/18 meningiomas and 5/50 melanomas, and, as expected, 2/2 giant cell tumors of the bone, and 7/7 xanthogranulomas (3).

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