**Monoclonal Mouse
Anti-Human
Epidermal Growth Factor Receptor**
Clone E30
Code M7239

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Epidermal Growth Factor Receptor, Clone E30, is intended for use in immunohistochemistry (IHC). Results aid in the classification of carcinomas expressing high levels of EGFR (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Synonyms for antigen	ErbB-1, HER1.
Summary and explanation	<p>Epidermal Growth Factor Receptor (EGFR) is a 170 kDa transmembrane glycoprotein belonging to a family of receptor tyrosine kinases comprising four members, EGFR (ErbB-1, HER1), ErbB-2 (HER2/neu), ErbB-3 (HER3), and ErbB-4 (HER4). These receptors are composed of an extracellular binding domain, a transmembrane lipophilic segment, and an intracellular protein tyrosine kinase domain (deficient in HER3), with a regulatory carboxyl-terminal segment. They become activated upon homo- or heterodimerization, which could be promoted either by ligand binding or by high receptor density due to overexpression. The ligands all contain an epidermal growth factor-like domain, and at least six different ligands have been identified for the EGFR. Activation by dimerization results in receptor autophosphorylation leading to recruitment of intracellular signalling pathways which are involved in the regulation of cell proliferation (2, 3). EGFR may be found expressed in a wide range of normal epithelial tissues and in tumors derived thereof (1, 4).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody supplied in liquid form as purified IgG from ascitic fluid. In 0.05 mol/L Tris-HCl, 15 mmol/L Na₂S₂O₃, 1% bovine serum albumin, pH 7.2.</p> <p>Clone: E30. Isotype: IgG1, kappa.</p> <p><u>Mouse IgG concentration</u>: see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Purified, denatured epidermal growth factor receptor.
Specificity	In CHO cells transfected with cDNA encoding wild type EGFR or EGFRvIII protein, membranes are strongly labeled by the antibody, whereas non-transfected CHO cells are negative.
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (Na₃N), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections</u>: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with proteinase K, trypsin (5) or pronase E (6) is required. Heat-induced epitope retrieval of tissues with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9 was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations</u>: The antibody can be used for labeling of acetone-fixed, frozen sections (5). The user must validate the staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution</u>: Monoclonal Mouse Anti-Human Epidermal Growth Factor Receptor (EGFR), Code M7239, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of tonsil or planocellular carcinoma and using 5 minutes proteolytic epitope retrieval with proteinase K (Code S3020), and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control</u>: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p>

Visualization: Dako EnVision+ /HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Staining interpretation Cells labeled by the antibody display staining predominantly of the cell membrane (5).

Performance characteristics **Normal tissues:** In normal oral mucosa, the antibody gives intense membranous labeling of epithelial cells. The labeling diminishes progressively from the basal towards the most superficial cell layer (5).

Abnormal tissues: The antibody labeled 21/111 cases of breast cancer (6). In glioblastoma, the antibody showed medium to strong labeling in about 75% of 55 cases (7). In 100 cases of oral squamous cell carcinomas divided in grade 1, 2 and 3 according to a patho-anatomical grading of the tumor front, the antibody only labeled the basal layer of the epithelium in grade 1, intensively labeled many tumor cells in grade 2 and nearly all tumor cells in grade 3 (8).

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
1. Gullick WJ. Prevalence of aberrant expression of the epidermal growth factor receptor in human cancers. Br Med Bull 1991;47:87-98.
2. Mendelsohn J, Baselga J. The EGF receptor family as targets for cancer therapy. Oncogene 2000;19:6550-65.
3. Schlessinger J, Schreiber AB, Levi A, Lax I, Libermann T, Yarden Y. Regulation of cell proliferation by epidermal growth factor. CRC Crit Rev Biochem 1983;14:93-111.
4. Gusterson B, Cowley G, Smith JA, Ozanne B. Cellular localisation of human epidermal growth factor receptor. Cell Biol Int Rep 1984;8:649-58.
5. Li T-J, Browne RM, Matthews JB. Expression of epidermal growth factor receptors by odontogenic jaw cysts. Virchows Arch A Pathol Anat Histopathol 1993;423:137-44.
6. Schroeder W, Biesterfeld S, Zillessen S, Rath W. Epidermal growth factor receptor – Immunohistochemical detection and clinical significance for treatment of primary breast cancer. Anticancer Res 1997;17:2799-802.
7. Schober R, Bilzer T, Waha A, Reifenberger G, Wechsler W, von Deimling A, et al. The epidermal growth factor receptor in glioblastoma: genomic amplification, protein expression, and patient survival data in a therapeutic trial. Clin Neuropathol 1995;14:169-74.
8. Störkel S, Reichert T, Reiffen KA, Wagner W. EGFR and PCNA expression in oral squamous cell carcinomas - a valuable tool in estimating the patient's prognosis. Oral Oncol, Eur J Cancer 1993;29B:273-7.

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