**Monoclonal Mouse
Anti-Human p53 Protein**
Clone DO-7

Code M7001

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human p53 Protein is intended for use in immunohistochemistry (IHC). The antibody labels wild-type and mutant-type p53 protein and is a useful tool for investigation of p53 accumulation in human neoplasias (1, 2). Differential classification of tumors is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>p53 is a nuclear phosphoprotein with a molecular mass of 53 kDa. Wild-type p53 protein is present in a wide variety of normal cells, but the protein has a very short half-life and thus is present in only minute amounts (1), generally below the detection level of immunocytochemical methods (3). Somatic mutation of the <i>p53</i> gene is a very frequent event in the development of human neoplasia, and because mutant p53 proteins often are much more stable than wild-type p53 protein, the mutant p53 protein accumulates to a high level (1). As examples, p53 protein accumulation was observed in 76% of 212 human malignant lesions, including breast, colon and stomach carcinomas, melanoma, embryonal carcinoma of the testis, transitional carcinoma of the urinary bladder, uterine carcinoma and soft tissue sarcomas (4).</p> <p>Wild-type p53 protein functions as a transcription factor, i.e. as a modulator which can turn crucial genes either on or off. It also inhibits DNA replication and is a check-point control molecule for progression of the cell cycle. Furthermore, p53 protein is involved in the regulation of apoptosis. In transfection assays, wild-type <i>p53</i> behaves as a tumor suppressor, while mutant <i>p53</i> behaves as a dominant transforming oncogene (1).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> DO-7 (1). <u>Isotype:</u> IgG2b, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Recombinant human wild-type p53 protein (1).
Specificity	<p>SDS-PAGE analysis of immunoprecipitates formed between lysate of the BT474 breast cancer cell line and the antibody shows reaction with a 53 kDa protein corresponding to p53 (1).</p> <p>In Western blotting of lysate of the A431 human vulval carcinoma cell line, the antibody labels a 53 kDa band, corresponding to the mutant type p53, which is expressed by A431. The epitope recognized by the antibody is located between the N-terminal amino acids 1 and 45 and possibly between amino acids 37 and 45 of the human p53 protein (1).</p> <p>In IHC the antibody labels mutant-type p53 in the A431 cell line and wild-type p53 in the SVK14 cell line (SV40-transformed keratinocyte line) (1).</p>
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin or methacarn (1). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed, frozen sections and acetone-fixed, cultured cells (1). The user must validate the staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human p53 Protein, Code M7001, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of breast carcinoma and epidermoid carcinoma cell line and using 20 minutes heat-induced</p>

SSM7001CEEFG_01 p. 1/4

epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0 and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG2b, Code X0944, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.

Visualization: Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.

Staining interpretation	Cells labeled by the antibody generally display a nuclear staining pattern, but cytoplasmic staining has been reported in some cases (5).
Performance characteristics	<p>Normal tissues: In normal and reactive mesothelium the antibody labeled 0/40 cases, and in 27 mesotheliomas, normal cells, e.g. fibroblasts and endothelial cells were not labeled (2).</p> <p>Abnormal tissues: In follicular lymphoma an increasing accumulation of p53 in centroblasts is observed with morphological progression resulting in 1/16 cases of grade I, 10/21 cases of grade II, and 6/6 cases of grade III being labeled (3). In mesotheliomas the antibody labeled 7/26 cases of epithelial type (1 to 25% labeled cells), 1/7 cases of mixed type (25 to 50% labeled cells), and 1/3 cases of mesenchymal type (more than 75% labeled cells) (2).</p>

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