**Monoclonal Mouse  
Anti-Human CD4  
Clone 4B12**

**Code M7310**

**ENGLISH**

**Intended use** For in vitro diagnostic use.  
Monoclonal Mouse Anti-Human CD4, Clone 4B12, is intended for use in immunohistochemistry (IHC). This antibody labels thymocytes and T-helper cells (1). Results aid in the classification of anaplastic large cell lymphoma (2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

**Summary and explanation** CD4 is a 55 kDa glycoprotein with five immunoglobulin-like external domains, a transmembrane domain and a highly conserved intracellular domain. CD4 acts as a co-receptor in MHC class II-restricted antigen-induced T-cell activation and may increase cell adhesion (1). The CD4 molecule also serves as a receptor for the human immunodeficiency virus on T cells and monocytes/macrophages (4).  
CD4 is expressed during T-cell development but is absent in immature thymocytes (3). It is expressed in T-helper cells, about 80-90% of mature thymocytes and 55-65% of mature peripheral T cells. It is also present on a subset of suppressor or cytotoxic T cells. CD4 is also expressed on monocytes/macrophages, Langerhans' cells and other dendritic cells, but it is not expressed on B cells (1).

CD4 plays an important role in immunophenotyping of reactive lymphocytes and lymphoproliferative disorders. The majority of peripheral T-cell lymphomas are derived from the helper T-cell subset (1). Antibodies to CD4 may also aid in the classification of mycosis fungoides and unspecified peripheral T-cell lymphoma (3-5).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

**Reagent provided** Monoclonal mouse antibody provided in liquid form as cell culture supernatant in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.  
Clone: 4B12. Isotype: IgG1, kappa.  
Mouse IgG1 concentration mg/L: See label on vial.  
The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

**Immunogen** Recombinant protein corresponding to the external domain of the CD4 molecule.

- Precautions**
- For in vitro diagnostic use.
  - For professional users.
  - This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
  - As with any product derived from biological sources, proper handling procedures should be used.
  - Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
  - Unused solution should be disposed of according to local, State and Federal regulations.

**Storage** Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

**Specimen preparation** Paraffin sections:  
The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin.  
Pre-treatment of deparaffinized\* tissue sections with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues with HIER using Dako EnVision FLEX Target Retrieval Solution, High pH (50x) (Code K8004).  
Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link. For details, please refer to PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove slide rack from PT tank and immediately dip slides in jar/tank (e.g., PT Link Rinse Station, Code PT109) containing diluted room temperature EnVision FLEX Wash Buffer (20X) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes.  
\*Paraffin-embedded sections: As alternative specimen preparation, deparaffinization, rehydration and epitope retrieval can be performed in the PT Link using a modified procedure. See the PT Link User Guide for instructions. After the staining procedure has been completed, the sections must be dehydrated, cleared and mounted using a permanent mounting method.  
The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.

**Staining procedure** These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.  
Dilution: Monoclonal Mouse Anti-Human CD4, Code M7310, may be used at a dilution range of 1:40-1:80 when applied on pretreated, formalin-fixed, paraffin-embedded sections using a 20 minute incubation at room temperature. It is recommended that the antibody is diluted in EnVision FLEX Antibody Diluent (Codes K8006) or Dako Antibody Diluent (Code S0809). The recommended negative control reagent is Dako Negative Control, Mouse IgG1 (Code X0931), diluted to the same IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately prior to use.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include tonsils and liver and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics".

Visualization: EnVision FLEX, High pH (Code K8000/K8010), using a 20 minute incubation at room temperature is recommended. Follow the procedure enclosed with the selected visualization system(s).

**Staining interpretation** Cells labeled by the antibody display membrane staining.

**Performance characteristics** Normal tissues: In tonsil, the crowded and isolated T-helper cells show a moderate to strong staining reaction. In liver, the Kupffer and endothelial cells of the sinusoids show a weak to moderate staining reaction.

Abnormal tissues: The antibody labeled 1/1 mucosis fungoides, and 2/2 peripheral T-cell lymphomas. In a study of 75 pediatric patients, the antibody labeled 22/51 anaplastic large cell lymphomas (2).

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- Dabbs DJ. Ed. Diagnostic Immunohistochemistry, second ed 2006 Elsevier Inc; p. 141-50.