**Negative Control
Mouse IgG1
Code X0931**

ENGLISH	
Intended use	For in vitro diagnostic use. X0931 is intended for use in flow cytometry and immunocytochemistry. X0931 has been produced especially as a negative control for Dako monoclonal mouse antibodies of isotype IgG1 and provided as cell culture supernatants.
Summary and explanation	Mouse isotype-specific negative control reagents are typically monoclonal antibodies directed against non-human antigens. These reagents measure the amount of primary antibody binding that is unrelated to the target antigen of the human cell. This non-immunological binding may be caused by Fc receptors on the cells. The isotypic control should be of the same isotype as the monoclonal mouse antibody used as the test reagent.
Reagent provided	X0931 is a monoclonal mouse antibody, clone DAK-GO1, provided in liquid form as cell culture supernatant (RPMI 1640 medium, containing calf serum), dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG1 concentration mg/L:</u> See label on vial.
Specificity	The antibody produced by clone DAK-GO1 is directed towards <i>Aspergillus niger</i> glucose oxidase, an enzyme which is neither present nor inducible in mammalian tissues. The antibody does stain prostate epithelial cells at high concentrations and as such should not be used at concentrations higher than 1:100 in this tissue. The antibody does stain pancreas cells at high concentrations and as such should not be used at concentrations higher than 1:200 in this tissue.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Dilution guidelines	IMMUNOCYTOCHEMISTRY X0931, should be adjusted to the same mouse IgG1 concentration as the test antibody. Unless the stability of the diluted negative control has been established in the actual staining procedure, it is recommended to dilute the reagent in Dako background reducing Antibody Diluent, Code S3022.
Dilution guidelines	FLOWCYTOMETRY X0931 may be used at a volume of 10 µL when adjusted to the same mouse IgG1 concentration as the test antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory.

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