**Polyclonal Rabbit
Anti-Goat Immunoglobulins/Biotinylated**
Code No./ Code/ Code-Nr. E 0466
Edition/ Ausgabe 14.11.02

ENGLISH

Intended use	For in vitro diagnostic use. E 0466 is intended for use in immunocytochemistry, ELISA and immunoblotting. Owing to a particularly strong cross-reaction with sheep immunoglobulins, E 0466 may serve as an excellent link antibody between goat or sheep primary antibodies and avidin or streptavidin labeled with enzymes or other markers. Cross-reaction with human immunoglobulins is negligible.
Reagent provided	Biotinylated, affinity-isolated rabbit anti-goat immunoglobulins provided in liquid form in 0.05 mol/L Tris/HCl, 0.1 mol/L NaCl, 15 mmol/L NaN ₃ , pH 7.2. <u>Concentration of specific antibody g/L:</u> see label on vial.
Preparation	<ol style="list-style-type: none">1. The antibody has been solid-phase absorbed with human serum proteins.2. The absorbed antibody has been further purified by affinity isolation on a column with insolubilized goat immunoglobulins using an elution procedure which ascertains a good recovery of high-affinity antibodies.3. Finally, biotinylation has been performed with activated biotin containing an aminocaproic acid spacer arm.
Immunogen	Immunoglobulins, mainly IgG, isolated from goat serum.
Specificity	The antibody reacts with goat immunoglobulins of all classes. Cross-reactions as determined by ELISA: Human immunoglobulins, less than 1%. Fetal calf serum, less than 2%. Mouse and rat immunoglobulins, about 25%. Ox immunoglobulins, about 50%. Sheep immunoglobulins, well over 80% (goat and sheep immunoglobulins have many common epitopes).
Precautions	<ol style="list-style-type: none">1. For professional users.2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. If unexpected reactions are observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact our Technical Services.
Dilution guidelines	Immunocytochemistry: 1:400-1:800, when used as link antibody between a goat or sheep primary antibody and avidin or streptavidin labelled with enzyme or fluorochrome. Immunoblotting (25 ng of goat or sheep immunoglobulins): 1:5 000-1:20 000. ELISA: 1:10 000-1:50 000. Unless the stability in the actual test system has been established, it is recommended to dilute the reagent immediately before use.

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