**Monoclonal Mouse  
Anti-Human  
Cytokeratin 10**  
Clone DE-K10  
Code M7002

**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use.  Monoclonal Mouse Anti-Human Cytokeratin 10, Clone DE-K10, is intended for use in immunohistochemistry (IHC). The antibody labels all suprabasal cells in the epidermis and is a useful aid for the classification of more differentiated squamous cell carcinomas (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.
<b>Summary and explanation</b>	The cytokeratins (CKs) belong to the intermediate filaments, which create a cytoskeleton in almost all eukaryotic cells. In contrast to other intermediate filaments, CKs are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. CKs are generally held to belong to the most fundamental markers of epithelial differentiation, and until now, 20 distinct CK polypeptides have been revealed in various human epithelia (3). The CKs can be divided into an acidic type A (class I) and a neutral-basic type B (class II) subfamily. CK 10 is an intermediate-sized, acidic type I cytokeratin, with a molecular mass of 56.5 kDa, expressed in all suprabasal cells simultaneously with CK 1 but is absent in basal cells. Together, CK 1 and CK 10 represent some of the first markers of epidermal differentiation (1, 4). CK 10 has been implicated as marker for more differentiated cells in squamous cell carcinomas (1, 2).  Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> . <u>Clone:</u> DE-K10 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.  The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
<b>Immunogen</b>	Cytoskeletal preparation extracted from human ectocervical epithelium (1).
<b>Specificity</b>	In Western blotting of cytoskeletal preparations from human epidermis, the antibody labels a single 56.5 kDa band corresponding to CK 10 (1). In two-dimensional immunoblotting of cytoskeletal proteins from human skin, the antibody labels a single dot corresponding to CK 10 (1).
<b>Precautions</b>	<ol style="list-style-type: none"><li>For in vitro diagnostic use.</li><li>For professional users.</li><li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li><li>As with any product derived from biological sources, proper handling procedures should be used.</li><li>Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li><li>Unused solution should be disposed of according to local, State and Federal regulations.</li></ol>
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is not required, but optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0 or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. <u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. <u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed, frozen sections (1). The user must validate the staining procedure.
<b>Staining procedure</b>	These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms. <u>Dilution:</u> Monoclonal Mouse Anti-Human Cytokeratin 10, Code M7002, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human epidermis and using 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. <u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

<b>Product-specific limitations</b>	The antibody labels all the cells in the sebaceous glands in frozen sections, but only some cells in the basal layer of the acini on formalin-fixed, paraffin-embedded sections (1).
<b>Staining interpretation</b>	Cells labeled by the antibody display a cytoplasmic staining pattern (1).
<b>Performance characteristics</b>	<u>Normal tissues:</u> The antibody labels all suprabasal layers of the epidermis, thymic Hassall's bodies and a variable amount of cells in some non-cornyfing stratified epithelia including vagina, ectocervix and tongue (1). Some cells are also labeled in hair follicles, sebaceous glands, and in the inner layer of the large secretory ducts of sweat glands (1). All simple and glandular epithelia and transitional epithelium of urinary bladder and urethra are not labeled by the antibody (1). <u>Abnormal tissues:</u> The antibody labeled 21/26 (81%) squamous cell carcinomas of the vulva, mostly in more differentiated parts (1). Various adenocarcinomas were not labeled by the antibody (1). Among 268 cervical biopsies, including 55 normal, 20 low grade squamous intraepithelial lesions (LG-SIL), 8 high grade squamous intraepithelial lesions (HG-SIL) and 29 invasive carcinomas, the median percentage of immunostaining for CK 10 was 43% in normal tissue, 41% in LG-SIL, 17% in HG-SIL and 1% in invasive carcinomas (2).

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


<b>References/ Références/ Literatur</b>												
<ol style="list-style-type: none"><li>Ivanyi D, Asink A, Groeneveld E, Hageman PC, Mooi WJ, Heintz APM. New monoclonal antibodies recognizing epidermal differentiation-associated keratins in formalin-fixed, paraffin-embedded tissue. Keratin 10 expression in carcinoma of the vulva. J Pathol 1989;159:7-12.</li> <li>Maddox P, Sasieni iiP, Szarewski A, Anderson M, Hanby A. Differential expression of keratins 10, 17, and 19 in normal cervical epithelium, cervical intraepithelial neoplasia, and cervical carcinoma. J Clin Pathol 1999;52:41-6.</li> <li>Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Hermann, Harris, editors. Subcellular Biochemistry. Volume 31, New York: Plenum Press;</li> <li>Moll R, Franke WW, Schiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. Cell 1982;31:11-24.</li> <li>Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. Subcellular Biochem; Intermediate Filaments. Herrmann, Harris, editors. Plenum Press, New York; 1998; Volume 31. p. 205-60.</li></ol>												
<b>Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole</b>												
<table border="1"> <tbody><tr> <td><b>REF</b></td> <td>Catalogue number Référence du catalogue Bestellnummer</td> <td> <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p></td> <td> <p>Manufacturer Fabricant Hersteller</p></td> </tr> <tr> <td><b>IVD</b></td> <td>In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum</td> <td><b>LOT</b></td> <td>Batch code Code du Lot Chargenbezeichnung</td> </tr> <tr> <td></td> <td>Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten</td> <td></td> <td>Use by Utiliser jusque Verwendbar bis</td> </tr> </tbody></table>	<b>REF</b>	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>	<b>IVD</b>	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	<b>LOT</b>	Batch code Code du Lot Chargenbezeichnung		Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis
<b>REF</b>	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>									
<b>IVD</b>	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	<b>LOT</b>	Batch code Code du Lot Chargenbezeichnung									
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis									

Revision 2017.02
<p> </p>