**Monoclonal Mouse
Anti-Human
Cytokeratin 7
Clone OV-TL 12/30
Code M7018**

ENGLISH**Intended use**

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human Cytokeratin 7, Clone OV-TL 12/30, is intended for use in immunohistochemistry (IHC). The antibody labels glandular and transitional epithelial cells and is a useful aid for the classification of adenocarcinoma of the lung (1, 2), breast and endometrium, thyroid gland (1, 3) and ovary (1, 4), as well as transitional cell (urothelial) carcinoma (1), and chromophobe renal cell carcinoma (5). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

Summary and explanation

Cytokeratin 7 belongs to the intermediate filaments, which create a cytoskeleton in almost all eukaryotic cells. In contrast to other intermediate filaments, cyto keratins (CKs) are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. CKs are generally held to belong to the most fundamental markers of epithelial differentiation, and until now, 20 distinct CK polypeptides have been revealed in various human epithelia (6, 7). The CKs can be divided into an acidic type A (class I) and a neutral-basic type B (class II) subfamily. CK7, a 54 kDa protein, belongs to the neutral-basic type B subfamily, and its distribution is confined to glandular and transitional epithelia (6).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required; Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

Reagent provided

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCL, pH 7.2, and containing 15 mmol/L NaN₃.

Clone: OV-TL 12/30 (1). Isotype: IgG1, kappa.

Mouse IgG concentration: see label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

OTN 11, ovarian carcinoma cell line (1, 8).

Specificity

In one- and two-dimensional SDS-PAGE and immunoblotting of Triton X-100-extracted cytoskeleton preparations of several cell cultures and tissues, the antibody labeled a 54 kDa band corresponding to cytokeratin 7 (1).

Precautions

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin sections: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin or Bouin's (3). Pre-treatment of deparaffinized tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval of tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labeling frozen sections (1, 3) and fixed cell smears (9). The user must validate the staining procedure.

Staining procedure

These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.

Dilution: Monoclonal Mouse Anti-Human Cytokeratin 7, Code M7018, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human breast and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0 and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen.

Visualization: Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Product-specific limitations

Exceptions to the generally expected reactivity pattern may occur, for example, CK 7-labeled hepatocytes have been observed in patients with acute and chronic cholestasis (1).

Cells labeled by the antibody display a cytoplasmic staining pattern.

Performance characteristics

Normal tissues: The antibody consistently labels a large number of simple-, complex- and transitional epithelia, including biliary and pancreatic ducts, lung alveoli, endometrium, distal convoluted tubules and collecting ducts of the kidney (simple epithelia), bronchial and bronchiolar epithelium, ducts of prostate, luminal cells of fallopian tube and endocervix, bronchial-, breast-, salivary-, sweat- and endocervical glands, placental trophoblasts (complex epithelia), and all cell layers of urothelium (transitional epithelium) (3). Additionally, the antibody labels ovarian mesothelium (4), and labeling has also been observed in luminal- and basal cells of the prostate and myoepithelial cells (3). Further, in frozen sections, the antibody has been shown to label the rete epithelium in the testis, epididymis epithelium, and the surface epithelium of the stomach and duodenum (1). Non-epithelial tissues, such as connective tissue, blood vessels, and lymphoid tissue, are negative with the antibody (1).

Abnormal tissues: In the human ovary, 12/12 cystoma simplex, 12/12 cystadenomas and 60/60 carcinomas were labeled by the antibody (4). Furthermore, 6/6 endometrial, 4/4 endocervical-, 3/3 breast-, and 3/3 thyroid carcinomas were labeled (3). In frozen tissues of the female genital tract, 1/1 Brenner's ovarian tumor, 21/21 serous- and 6/6 mucinous cystadenocarcinomas, 4/4 nonclassified ovarian carcinomas, 8/8 endometrioid carcinomas, 1/1 adenosquamous carcinoma of the cervix, 1/1 metastasis of fallopian tube adenocarcinoma in the cervix, and 1/1 placental choriocarcinoma revealed labeling with the antibody (1). In the human lung, 20/20 adenocarcinomas of different grades, including 4 bronchioalveolar carcinomas, were labeled with the antibody, whereas 24/24 cases of squamous cell carcinomas were unlabeled, as were 6/6 cases of large cell anaplastic carcinomas (2) and 10/10 cases of mesotheliomas (9). In the human gastrointestinal tract, 5/6 differently graded adenocarcinomas and 3/3 poorly differentiated signet ring cell carcinomas of the stomach, and 1/1 carcinoma of the pancreas were labeled by the antibody, whereas anaplastic stomach-, small intestine-, and rectal carcinomas as well as colon adenocarcinomas of different grades were unlabeled (3). In frozen tissues, 1/1 anaplastic adenocarcinoma of the stomach, 1/1 adenocarcinoma of jejunum, 2/2 intestinal carcinoids, 1/1 hepatoblastoma and 3/8 hepatocellular carcinomas were labeled. Also 4/4 transitional cell carcinomas were labeled (1). 6/6 chromophobe renal cell carcinomas were labeled, whereas 8/11 oncocytomas were unlabeled with the antibody (5).

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

- Ramaekers F, van Niekerk C, Poels L, Schaafsma E, Huijsmans A, Robben H, et al. Use of monoclonal antibodies to keratin 7 in the differential diagnosis of adenocarcinomas. *Am J Pathol* 1990;136:641-55.
- van de Molengraft FJMM, van Niekerk CC, Jap PHK, Poels LG. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. *Histopathology* 1993;22:35-8.
- van Niekerk CC, Jap PHK, Ramaekers FCS, van de Molengraft F, Poels LG. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. *J Pathol* 1991;165:145-52.
- van Niekerk CC, Boerman OC, Ramaekers FCS, Poels LG. Marker profile of different phases in the transition of normal human ovarian epithelium to ovarian carcinomas. *Am J Pathol* 1991;138:455-63.
- Leroy X, Moukassa D, Copin MC, Saint F, Mazeman E, Gosselin B. Utility of cytokeratin 7 for distinguishing chromophobe renal cell carcinoma from renal oncocytoma. *Eur Urol* 2000;37:484-7.
- Moll R, Franke WW, Schiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: Patterns of expression in normal epithelia, tumors and cultured cells. *Cell* 1982;31:11-24.
- Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Hermann, Harris, editors. *Subcellular biochemistry*. Volume 31, New York: Plenum Press; 1998. p. 205-62.
- Poels LG, Jap PHK, Ramaekers FCS, Scheres JMJC, Thomas CMG, Vooijs PG, et al. Characterization of a hormone-producing ovarian carcinoma cell line. *Gynecol Oncol* 1989;32:203-14.
- Baars JH, de Ruijter JLM, Smedts F, van Niekerk CC, Poels LG, Seldenrijk CA, et al. The applicability of a keratin 7 monoclonal antibody in routinely Papanicolaou-stained cytologic specimens for the differential diagnosis of carcinomas. *Am J Clin Pathol* 1994;101:257-61.

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