**Monoclonal Mouse
Anti-Human
Cytokeratin
Clone MNF116
Code M0821**

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Cytokeratin, Clone MNF116, is intended for use in immunohistochemistry (IHC). The antibody labels epithelial tissues from simple glandular to stratified squamous epithelium, and is a useful aid for the classification of neoplasms of epithelial origin (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>The cytokeratins (CKs) belong to the intermediate filaments, which create a cytoskeleton in almost all eukaryotic cells. In contrast to other intermediate filaments, CKs are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. CKs are generally held to belong to the most fundamental markers of epithelial differentiation, and 20 distinct CK polypeptides have been revealed in various human epithelia (3). The CKs can be divided into an acidic type A (class I) and a neutral-basic type B (class II) subfamily.</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> MNF 116. <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Crude extract of splenic cells from a nude mouse engrafted with MCF-7 cells (human breast carcinoma cell line).
Specificity	The antibody is a broad spectrum anti-keratin reagent reacting with intermediate and low-molecular-weight keratins. Thus in immunoblotting, the antibody labels a number of discrete bands ranging from 40 to 58 kDa, corresponding to cytokeratins 5, 6, 8, 17 and probably 19.
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with trypsin (1, 4), pepsin (2), proteinase K or heat-induced epitope retrieval is required. For heat-induced epitope retrieval, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human Cytokeratin, Code M0821, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p> <p><u>Visualization:</u> Dako EnVision+/-HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>
Staining interpretation	Cells labeled by the antibody display a cytoplasmic staining pattern.

Performance characteristics Normal tissues: The antibody labels a number of different epithelia, e.g. simple epithelia in biliary and pancreatic ducts, proximal and distal convoluted renal tubules, thyroid-, parathyroid-, gastric-, gall bladder-, small intestinal-, large intestinal-, and endometrial epithelium, as well as hepatocytes (1). The epithelial lining of the rete ovarii, oophoron, Fallopian tube and ovarian surface is also reactive with the antibody (4). In stratified squamous epithelium of skin, the antibody labels the basal layer, as well as the outer root sheath of follicular epithelium, and the basal layer of sebaceous, eccrine, and apocrine glands and ducts (2). Non-squamous stratified epithelia of esophagus and ectocervix are also reported positive with the MNF 116 antibody (1). Transitional epithelium (urothelium) and different complex epithelia, e.g. ducts and acini of breast, endocervical glands, prostate epithelium, epididymis, bronchial epithelium as well as ducts and acini of salivary glands are also labeled by the antibody (1). Mesenchymal cells are generally not reactive with the antibody, however a weak and focal labeling of smooth muscle cells is reported. Endothelium, skeletal muscle, cartilage, and lymphoid tissues are not labeled by the antibody (1).

Abnormal tissues: Of epithelial tumors, the antibody labeled 4/4 colorectal carcinomas, 3/3 gastric carcinomas, 4/4 breast carcinomas, 3/3 prostatic carcinomas, 3/3 renal cell carcinomas, 3/4 hepatocellular carcinomas, 3/3 transitional cell carcinomas, 2/3 carcinoids of appendix, and, moreover, 2/2 teratomas and 2/2 pleomorphic adenomas (epithelial elements), and 14/14 squamous cell carcinomas, i.e. 6 of epidermal-, 3 of cervical- and 5 of bronchial origin (1). In addition, 16/16 cases of mixed tumors and myoepitheliomas arising in soft tissues, displayed immunoreactivity with the antibody (5). Labeling of mesothelioma has also been observed. The antibody, further, labels micrometastatic cells of non-small cell lung carcinomas in lymph nodes (6). Non-epithelial tumors, e.g. melanomas, lymphomas, benign and malignant fibrous histiocytomas, leiomyomas, liposarcomas, chondrosarcomas, Ewings sarcomas, angiofibromas, angiosarcomas, atypical fibroxanthomas, juvenile xanthogranulomas, hemangiomas, granular cell tumors, hair follicle myxomas, epithelioid histiocytomas and dermatofibrosarcomas are not labeled by the antibody. However, in 1/6 cases (1), and 1/2 cases (2) of leiomyosarcomas, the antibody displayed a weak and focal reactivity.

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
- Goddard MJ, Wilson B, Grant JW. Comparison of commercially available cytokeratin antibodies in normal and neoplastic adult epithelial and non-epithelial tissues. J Clin Pathol 1991;44:660-3.
- Prieto VG, Lugo J, McNutt NS. Intermediate- and low-molecular-weight keratin detection with the monoclonal antibody MNF116. An immunohistochemical study on 232 paraffin-embedded cutaneous lesions. J Cutan Pathol 1996;23:234-41.
- Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Hermann, Harris, editors. Subcellular Biochemistry. Volume 31, New York: Plenum Press; 1998. p. 205-62.
- Woolnough E, Russo L, Khan MS, Heatley MK. An immunohistochemical study of the rete ovarii and epoophoron. Pathology 2000;32:77-83.
- Kilpatrick SE, Hitchcock MG, Kraus MD, Calonje E, Fletcher CDM. Mixed tumors and myoepitheliomas of soft tissue. A clinicopathologic study of 19 cases with a unifying concept. Am J Surg Pathol 1997;21:13-22.
- Nicholson AG, Graham ANJ, Pezzella F, Agneta G, Goldstraw P, Pastorino U. Does the use of immunohistochemistry to identify micrometastases provide useful information in the staging of node-negative non-small cell lung carcinomas? Lung Cancer 1997;18:231-40.

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