**Monoclonal Mouse
Anti-Human
Smooth Muscle Actin
Clone 1A4
Code M0851**

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Smooth Muscle Actin, Clone 1A4, is intended for use in immunohistochemistry (IHC). The antibody labels smooth muscle cells, myofibroblasts and myoepithelial cells, and is a useful aid for classification of leiomyomas, leiomyosarcomas (1, 2), and pleomorphic adenomas (3). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>Cytoplasmic actins, which belong to the microfilament system of cytoskeleton proteins, are some of the most conserved eukaryotic proteins. Six isoforms of the actin protein exist. All isoforms possess very similar amino acid sequences, with no isoform sharing less than 93% identity with any other isoform. The N-terminal region appears to be least identical (4). Four isoforms, α-skeletal actin, α-cardiac actin, α-smooth actin, and γ-smooth actin, are mainly expressed in skeletal, cardiac, and smooth muscle while the two isoforms, β-cyto actin and γ-cyto actin are ubiquitously expressed (4).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> 1A4. The 1A4 clone is identical to the anti-asm-1 described in (5). <u>Isotype:</u> IgG2a, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	<p>N-terminal synthetic decapeptide of α-smooth muscle actin (5).</p>
Specificity	<p>In SDS-PAGE immunoblotting of the α-smooth muscle isoform of actin, the antibody labels a band corresponding to α-smooth muscle actin (5).</p>
Precautions	<ol style="list-style-type: none">1. For in vitro diagnostic use.2. For professional users.3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.4. As with any product derived from biological sources, proper handling procedures should be used.5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. However, Dako Target Retrieval Solution, Code S1700 was found inefficient. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed, frozen sections (1). The user must validate the staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human Smooth Muscle Actin, code M0851, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of normal human colon and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG2a, Code X0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p> <p><u>Visualization:</u> Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>
Staining interpretation	<p>Cells labeled by the antibody display a cytoplasmic staining pattern.</p>

Performance characteristics Normal tissues: The antibody labels smooth muscle cells in blood vessels and, additionally, salivary ducts and myoepithelial cells around acini in salivary glands (3). Smooth muscle cells in 35/36 normal uterine myometria were also positively labeled (2). Further, a temporal labeling of perisinusoidal liver cells has been observed (6). In frozen tissues, the antibody labels myofibroblasts and myoepithelial cells around acini and ducts of the breast, whereas epithelia (adeno, squamous), lymphocytes, cardiac- and skeletal muscle cells, endothelial cells, fat cells, Schwann cells and fibroblast are negative (1).

Abnormal tissues: The antibody labeled 24/26 leiomyomas, 6/7 atypical leiomyomas and 21/25 leiomyosarcomas of the uterus, as well as 13/13 extruterine nongastrointestinal spindle leiomyosarcomas (2). Moreover, the antibody labeled a variable amount of cells in 8/8 pseudosarcomatous myofibroblastic tumors of the urinary bladder in children (7). In pleomorphic adenomas, the antibody labeled tumor epithelial cells (myoepithelial cells) in 19/20 cases (3). In frozen tissues, the antibody, in addition to the labeling of 5/5 leiomyomas and 6/7 leiomyosarcomas, also labeled 4/22 malignant fibrous histiocytomas and 1/2 rhabdomyosarcomas. 6/6 malignant schwannomas were non-labeled, as also 13/13 other soft tissue tumors, including 1 fibrosarcoma, 6 liposarcomas, 1 angio-sarcoma, 1 capillary hemangioma, 1 Triton tumor, and 3 synovial sarcomas (1).

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

References / Bibliographie / Literaturnachweise

1. Roholl PJM, Elbers HRJ, Prinsen I, Claessens JAJ, van Unnik JAM. Distribution of actin isoforms in sarcomas: an immunohistochemical study. Hum Pathol 1990;21:1269-74.
2. Rizeq MN, van de Rijn M, Hendrickson MR, Rouse RV. A comparative immunohistochemical study of uterine smooth muscle neoplasms with emphasis on the epithelioid variant. Hum Pathol 1994;25:671-7.
3. Brennan PA, Umar T, Zaki GA, Langdon JD, Spedding A, Buckley J, et al. Are myoepithelial cells responsible for the widespread expression of inducible nitric oxide synthase in pleomorphic adenoma? An immunohistochemical study. J Oral Pathol Med 2000;29:279-83.
4. Perrin BJ and Ervasti JM. The actin gene family: Function follows isoform. Cytoskel 2010; 67:630-34.
5. Skalli O, Ropraz P, Trzeciak A, Benzouana G, Gillesen D, Gabbiani G. A monoclonal antibody against α -smooth muscle actin: a new probe for smooth muscle differentiation. J Cell Biol 1986;103:2787-96.
6. Schmitt-Gräff A, Krüger S, Bocharad F, Gabbiani G, Denk H. Modulation of alpha smooth muscle actin and desmin expression in perisinusoidal cells of normal and diseased livers. Am J Pathol 1991;138:1233-42.
7. Hojo H, Newton WA, Hamoudi AB, Qualman SJ, Wakasa H, Suzuki S, et al. Pseudosarcomatous myofibroblastic tumor of the urinary bladder in children: a study of 11 cases with review of the literature. Am J Surg Pathol 1995;19:1224-36.

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