**Monoclonal Mouse
Anti-Vimentin**
Clone Vim 3B4
Code M7020

ENGLISH

| | |
|--------------------------------|---|
| Intended use | For in vitro diagnostic use. Monoclonal Mouse Anti-Vimentin, Clone Vim 3B4, is intended for use in immunohistochemistry (IHC). The antibody labels primarily cells of mesenchymal origin in normal and neoplastic tissues, and is a useful aid for classification of tumors of mesenchymal origin (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains. |
| Summary and explanation | Vimentin is a 55 kDa intermediate filament (IF) protein, which form part of the cytoskeleton of vertebrate cells. Among the IFs, vimentin belongs to class III and is expressed in the cytoplasm of cells of mesenchymal origin (2). Being the predominant IF protein in mesenchymal cells vimentin was initially believed to be a useful marker for classifying non-epithelial neoplasms. The co-expression of intermediate filaments, particularly vimentin and cytokeratin, has been demonstrated in a variety of normal cells/tissues and in neoplastic lesions, necessitating the use of a panel of antibodies in tumor classification (1). Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations. |
| Reagent provided | Monoclonal mouse antibody provided in liquid form purified from cell culture supernatant. In 0.05 mol/L Tris-HCl, 1% bovine serum albumin (BSA) and 15 mmol/L Na ₂ S ₂ O ₄ , pH 7.2. <u>Clone:</u> Vim 3B4. <u>Isotype:</u> IgG2a, kappa. <u>Mouse IgG concentration:</u> See label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot. |
| Immunogen | Vimentin isolated from bovine eye lens (4). |
| Specificity | SDS-PAGE analysis of immunoprecipitates formed between the antibody and ³⁵ S-labeled proteins from metabolically labeled human osteogenic sarcoma 4-998 cells shows reaction with a ~60 kDa polypeptide, corresponding to vimentin. In addition, a few bands of lower molecular weight were observed and may represent degradation products of vimentin (4). The antibody reacts with an epitope that has been localized to the coil 2 part of the vimentin rod domain (4). As demonstrated by immunocytochemistry on frozen as well as formalin-fixed tissues, the antibody cross-reacts with the vimentin-equivalent protein in man (5-8). |
| Precautions | <ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations. |
| Storage | Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support. |
| Specimen preparation | <u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with proteinase K or heat-induced epitope retrieval is required. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone fixed, frozen sections (3, 6, 7). The user must validate the staining procedure. |
| Staining procedure | These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms. <u>Dilution:</u> Monoclonal Mouse Anti-Vimentin, Code M7020, may be used at a dilution range of 1:100-1:200 when applied on formalin-fixed, paraffin-embedded sections of malignant melanoma and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG2a, Code X0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. <u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit. <u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen. |

SSM7020CEEF01_p. 1/4

| | |
|------------------------------------|---|
| Staining interpretation | Cells labeled by the antibody display staining confined to the cytoplasm. |
| Performance characteristics | <u>Normal tissues:</u> The antibody strongly labels fibrocytes, lipocytes, smooth muscle cells, vascular endothelial cells, peripheral nerve (Schwann) cells, macrophages (including Kupffer cells) and myoepithelial cells of sweat and salivary glands. With variable intensity and distribution, the antibody also labels follicular cells of the thyroid, adrenal cortex, renal distal tubuli and mesangial and endothelial cells of the renal glomerulus in addition to pancreatic acinar cells (1). <u>Abnormal tissues:</u> Of carcinomas the antibody labeled 15/147 adeno-, 5/13 cholangio-, 2/29 hepatocellular, 1/6 small cell, 1/7 squamous cell and 3/10 undifferentiated. The antibody also labeled 16/18 melanomas, 4/4 meningiomas, 3/8 mesotheliomas, and 17/20 sarcomas. Of other lesions, the antibody labeled 1/10 carcinoids, 1/2 neuroblastomas, 4/7 paragangliomas, 1/1 pleomorphic adenoma, 3/3 schwannomas, and 1/5 thymomas (1). |

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
References / Bibliographie / Literaturnachweise

- Azumi N, Battifora H. The distribution of vimentin and keratin in epithelial and nonepithelial neoplasms. A comprehensive immunohistochemical study on formalin- and alcohol-fixed tumors. Am J Clin Pathol 1987;88:286-96.
- Herrmann H, Aebi U. Intermediate filaments and their associates: multi-talented structural elements specifying cytoarchitecture and cytodynamics. Curr Opin Cell Biol 2000;12:79-90.
- Herrmann H, Fouquet B, Franke WW. Expression of intermediate filament proteins during development of Xenopus laevis. I. cDNA clones encoding different forms of vimentin. Development 1989;105:279-98.
- Stathopoulos E, Naeve GS, Taylor CR, Epstein AL. LN-6: a monoclonal antibody to vimentin expressed in non-hematopoietic mesenchymal cells and derived tumors and reactive in B5-fixed, paraffin- embedded tissues. J Histochem Cytochem 1989;37:1363-70.
- Bohn W, Wieggers W, Beuttenmuller M, Traub P. Species-specific recognition patterns of monoclonal antibodies directed against vimentin. Exp Cell Res 1992;201:1-7.
- Olah I, Kendall C, Glick B. Anti-vimentin monoclonal antibody recognizes a cell with dendritic appearance in the chicken's bursa of Fabricius. Anat Rec 1992;232:121-5.
- Heid HW, Moll I, Franke WW. Patterns of expression of trichocytic and epithelial cytokeratins in mammalian tissues. I. Human and bovine hair follicles. Differentiation 1988;37:137-57.

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