**Polyclonal Rabbit  
Anti-Human IgA**
**Code A0262**
**ENGLISH**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>Polyclonal Rabbit Anti-Human IgA, Code A0262, is intended for use in immunohistochemistry (IHC). This antibody labels IgA-expressing cells in formalin-fixed, paraffin-embedded tissue. Differential classification of tumors is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
<b>Summary and explanation</b>	<p>The human immunoglobulins basically consist of two identical heavy chains (Mr about 50 000) and two identical light chains (Mr about 20 000). The four chains are covalently linked together by disulphide bonds. The light chains are either kappa or lambda. The five immunoglobulins, IgA, IgD, IgE, IgG and IgM, differ regarding the heavy chains, and IgA and IgM also occur in polymeric forms. The heavy chain of IgA is named the alpha-chain. In serum about 80% of IgA is monomeric, while in secretions, such as saliva, intestinal and bronchial mucus, nasal secretions, sweat, and breast milk and colostrum, the predominant form is the dimeric secretory IgA (sIgA). Besides 4 light-chains and 4 alpha-chains dimeric sIgA also contains a J-chain and a secretory component. The latter protects sIgA against secretory proteases. The Mr of dimeric sIgA is 390 000 (1).</p> <p>The normal B-cell population is polyclonal and expresses a range of different immunoglobulin molecules. In contrast, a majority of B-cell neoplasias are characterized by the proliferation of monoclonal cells expressing only one type of light chain, whereas more than one heavy chain can be expressed by the same cell. The restricted expression of immunoglobulins by monoclonal B-cell lineage proliferations makes antibodies specific for immunoglobulin light- and heavy chains useful for the classification of these neoplasias (2).</p> <p>Antibodies to IgA have been shown to aid in the classification of B-cell neoplasias (3-5).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
<b>Reagent provided</b>	<p>Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>.</p> <p><u>Protein concentration:</u> See label on vial. <u>Antibody titre (SRI):</u> 2200 mg/L (7).</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
<b>Immunogen</b>	IgA isolated from a pool of normal human sera.
<b>Specificity</b>	<p>The antibody reacts with the alpha-chains of human IgA. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins.</p> <p>The specificity of the antibody has been ascertained as follows:</p> <p><u>Crossed immunoelectrophoresis:</u> Only the IgA immunoprecipitation arch appears when using 12.5 µL antibody per square cm gel area against 2 µL human plasma. Staining: Coomassie Brilliant Blue.</p> <p><u>ELISA:</u> No significant reaction is seen in indirect ELISA using human IgG and IgM as coating antigens. In double antibody sandwich ELISA no significant reaction is observed with human plasma stripped of IgA.</p>
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For in vitro diagnostic use.</li> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> <li>Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
<b>Specimen preparation</b>	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with proteinase K or heat-induced epitope retrieval is required. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections (8). The user must validate the staining procedure.</p>
<b>Staining procedure</b>	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Polyclonal Rabbit Anti-Human IgA, Code A0262, may be used at a dilution range of 1:100-1:200 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Visualization:</u> Dako EnVision+ /HRP kits, e.g. Code K4009, are recommended. Follow the procedure enclosed with the selected visualization kit.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p>

<b>Staining interpretation</b>	Cells labelled by the antibody display staining of the cytoplasm and/or the cell membrane.
<b>Performance characteristics</b>	<p><b>Normal tissues:</b> In frozen sections of reactive lymph nodes, the antibody labels many plasma cells and, occasionally, follicle centre cells and intercellular material in the follicles (8).</p> <p><b>Abnormal tissues:</b> In two studies investigating 33 (5) and 29 (3) paraffin-embedded tissues from cases of plasma cell neoplasia of various types, the labelling with the antibody showed a good correlation with the monoclonal immunoglobulin found in serum or urine. A classification of reactive hyperplasia and monoclonal B-cell neoplasia was obtained when the antibody was used in combination with a panel of other antibodies (3, 6).</p>

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

## References/ Bibliographie/ Literaturnachweise

1. Klein J, Hořejší V. Immunology. 2nd ed. Abingdon (UK): Blackwell Science Ltd; 1999. p. 226-7, 238-46.
2. Leong ASY, Cooper K, Leong FJWA. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 217-9.
3. Taylor CR, Burns J. The demonstration of plasma cells and other immunoglobulin-containing cells in formalin-fixed, paraffin-embedded tissues using peroxidase-labelled antibody. J Clin Pathol 1974;27:14-20.
4. Taylor CR, Mason DY. The immunohistological detection of intracellular immunoglobulin in formalin-paraffin sections from multiple myeloma and related conditions using the immunoperoxidase technique. Clin exp Immunol 1974;18:417-29.
5. Taylor CR, Russell R, Chandor S. An immunohistologic study of multiple myeloma and related conditions, using an immunoperoxidase method. Am J Clin Pathol 1978;70:612-22.54.
6. Beschorner R, Horny H-P, Petrusch UR, Kaiserling E. Frequent expression of haemopoietic and non-haemopoietic antigens by reactive plasma cells: an immunohistochemical study using formalin-fixed, paraffin-embedded tissue. Histol Histopathol 1999;14:805-12.
7. Becker W. Determination of antisera titres using the single radial immunodiffusion method. Immunochem 1969;6:539-46.
8. Harris NL, Poppema S, Data RE. Demonstration of immunoglobulin in malignant lymphomas. Use of an immunoperoxidase technic on frozen sections. Am J Clin Pathol 1982;78:14-21.

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