**Monoclonal Mouse
Anti-Human CD3**
Clone F7.2.38

Code M7254

ENGLISH

Intended use

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human CD3, Clone F7.2.38, is intended for use in immunohistochemistry (IHC). The antibody labels CD3 and is a useful aid for the identification of T cells. Results aid in the classification of T-cell neoplasms (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

Synonyms for antigen

T3, CD3 complex (2).

Summary and explanation

The CD3 complex is composed of six polypeptides with usually four different transmembrane CD3 chains, γ (gamma), δ (delta), ϵ (epsilon), and ζ (zeta). Three different dimers, $\gamma\epsilon$, $\delta\epsilon$, and $\zeta\zeta$, constitute the CD3 complex. The Mr of CD3 ϵ is 20 000 (2).

CD3 is first detectable in early thymocytes and its appearance probably represents one of the earliest signs of commitment to the T-cell lineage (3). In cortical thymocytes, the antigen is predominantly present as an intracyto-plasmic constituent. It appears subsequently, at the medullar thymocyte stage, on the T-cell surface in close association with the T-cell receptor (TCR). Whereas the TCR forms the ligand-binding part of the TCR/CD3 complex, the function of the CD3 molecule is that of signalling, making CD3 a highly specific marker for T cells. No other cells are known to express the CD3 molecule, although two monoclonal antibodies raised against CD3 have been found to react with Purkinje cells in the cerebellum (4).

The CD3 molecule may be present in the great majority of T-cell neoplasms, although occasional tumors are encountered in which the antigen is lost as part of the neoplastic process (5). The CD3 molecule may also be expressed in some cases of malignant histiocytosis (6) and Hodgkin's lymphoma (7).

CD3 ϵ , and not the whole CD3 molecule, has been detected in the cytoplasm of natural killer (NK) cells (8), and CD3 ϵ may be a marker of nasal T-cell lymphomas which are thought to be of NK cell origin (9).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage, Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

Reagent provided

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN₃.

Clone: F7.2.38 (1). Isotype: IgG1, kappa.

Mouse IgG concentration: see label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

Purified CD3 $\epsilon\gamma\delta$ /CD3 ω (1).

Specificity

In Western blotting of FOX B cells transfected with CD3 ϵ cDNA (T116.4.3) or Jurkat E6.1 T cells, the antibody labels a 20 kDa band under reducing conditions, corresponding to CD3 ϵ . No labeling was observed in non-transfected FOX B cells (1).

The antibody is most likely directed against the cytoplasmic region of the CD3 ϵ chain and works equally well in immunohistochemistry as the Dako Polyclonal Rabbit Anti-Human CD3, Code A0452 (1).

Precautions

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin sections: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, Duboscq-Brasil, Bouin's, or Bouin Holland (1). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0 and pre-treatment of tissues with proteinase K. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labeling frozen sections (1). The user must validate the staining procedure.

Staining procedure

These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.

Dilution: Monoclonal Mouse Anti-Human CD3, Code M7254, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the

diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S 0809.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen.

Visualization: Dako EnVision+HRP kits, E.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Staining interpretation

Cells labeled by the antibody display membranous and/or cytoplasmic staining.

Performance characteristics

Normal tissues: In thymus the antibody strongly labels lymphoid cells in the medulla and cortex. In tonsil and lymph node, interfollicular areas containing the T lymphocytes are strongly labeled, as also small reactive T cells in the B-cell follicles. With the exception of scattered T cells, all other normal human tissues studied were negative with the antibody (1).

Abnormal tissues: The majority of T-cell lymphomas, 41/52 cases, were labeled by the antibody. Of non-Hodgkin's lymphomas, the antibody labeled 10/10 precursor T-lymphoblastic lymphomas/leukemias, 11/11 peripheral T-cell lymphomas, 8/8 mycosis fungoides/Sezary syndrome, 2/2 angioimmunoblastic T-cell lymphomas (AILD), 1/1 intestinal T-cell lymphoma, 6/14 anaplastic large cell lymphomas (ALCL), 3/5 extranodal NK/T-cell lymphomas, nasal type, and 0/1 T-cell large cell granular lymphocytic leukemia. In cases of Hodgkin's lymphoma, only 1/6 nodular sclerosis, and 1/1 unclassified showed focal cytoplasmic staining, whereas none of 5 cases of lymphocytic predominance, and 9 cases of mixed cellularity were labeled. None of 37 cases of different B-cell lymphomas were labeled (1).

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

References / Bibliographie / Literaturnachweise

- Alibaud L, Llobera R, Al Saati T, March M, Delsol G, Rubin B. A new monoclonal anti-CD3 antibody reactive on paraffin sections. J Histochem Cytochem 2000;48:1609-16.
- Jin YJ, Koyasu S, Moingeon P, Steinbrich R, Tarr GE, Reinherz EL. A fraction of CD3 epsilon subunits exists as disulfide-linked dimers in both human and murine T lymphocytes. J Biol Chem 1990;265:15850-3.
- Campana D, Thompson JS, Amlot P, Brown S, Janossy G. The cytoplasmic expression of CD3 antigens in normal and malignant cells of the T lymphoid lineage. J Immunol 1987;138:648-55.
- Garson JA, Beverley PCL, Coakham HB, Harper EI. Monoclonal antibodies against human T lymphocytes label Purkinje neurones of many species. Nature 1982;298:375-7.
- Picker LJ, Weiss LM, Medeiros LJ, Wood GS, Warnke RA. Immunophenotypic criteria for the diagnosis of non-Hodgkin's lymphoma. Am J Pathol 1987;128:181-201.
- Wilson MS, Weiss LM, Gatter KC, Mason DY, Dorfman RF, Warnke RA. Malignant histiocytosis. A reassessment of cases previously reported in 1975 based on paraffin section immunophenotyping studies. Cancer 1990;66:530-6.
- Cibull ML, Stein H, Gatter KC, Mason DY. The expression of the CD3 antigen in Hodgkin's disease. Histopathology 1989;15:597-605.
- Lanier LL, Chang C, Spits H, Phillips JH. Expression of cytoplasmic CD3ε proteins in activated human adult natural killer (NK) cells and CD3γ, δ, ε complexes in fetal NK cells. J Immunol 1992;149:1876-80.
- Ohno T, Yamaguchi M, Oka K, Miwa H, Kita K, Shirakawa S. Frequent expression of CD3ε in CD3 (Leu 4)-negative nasal T-cell lymphomas. Leukemia 1995;9:44-52.

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