**Monoclonal Mouse
Anti-Human
Ki-67 Antigen
Clone MIB-1**

Code M7240

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Ki-67 Antigen, Clone MIB-1, is intended for use in immunohistochemistry (IHC). The antibody is useful for the identification of the Ki-67 antigen in normal and neoplastic cells (1). Differential classification of tumors is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>The Ki-67 antigen is a nuclear protein, which is defined by its reactivity with monoclonal antibody from the Ki-67 clone (2). Two isoforms of 345 and 395 kDa have been identified (3). The Ki-67 antigen is preferentially expressed during all active phases of the cell cycle (G₁, S, G₂ and M-phases), but it is absent in resting cells (G₀-phase) (2). During interphase, the antigen can be exclusively detected within the nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes. The antigen is rapidly degraded as the cell enters the non-proliferative state (4), and there appears to be no expression of Ki-67 during DNA repair processes (5).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage, Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, 1% bovine serum albumin, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> MIB-1 (6). <u>Isotype:</u> IgG1, kappa. <u>Mouse Iq concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Human recombinant peptide corresponding to a 1002 bp Ki-67 cDNA fragment (6).
Specificity	In Western blotting of lysates of the multiple myeloma cell line, IM-9, the MIB-1 antibody labels bands of 345 and 395 kDa, identical to the bands labelled by the original Ki-67 antibody. Furthermore, Western blotting and competitive binding experiments clearly demonstrate that MIB-1, like the original Ki-67 antibody, reacts with an epitope encoded by a 66 bp repetitive element in the Ki-67 gene. In immunohistochemistry, the MIB-1 and the Ki-67 antibodies provide identical staining patterns on serial tonsillar frozen sections (6). The MIB-1 antibody recognizes native Ki-67 antigen and recombinant fragments of the Ki-67 molecule (6).
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin.</p> <p><u>Pre-treatment:</u> Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating deparaffinized tissues with HIER using diluted Dako Target Retrieval Solution, Low pH (10x) (Code S1699) or diluted EnVision FLEX Target Retrieval Solution, Low pH (50x) (Code K8005) for 20 minutes. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human Ki-67 Antigen, Code M7240, may be used at a dilution range of 1:75-1:150 when applied on formalin-fixed, paraffin-embedded sections of human tonsil or human intestinal mucosa and using 20 minutes heat-induced epitope retrieval in Target Retrieval Solution, Low pH (Code S1699/K8005), and 20 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality Control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p> <p><u>Visualization:</u> The recommended visualization system is EnVision FLEX, High pH (Code K8000/K8010) using a 20 minutes incubation at room temperature. Follow the procedure enclosed with the visualization kit.</p> <p>Note: Use diluted Dako Target Retrieval Solution, Low pH (10x) (Code S1699) or diluted EnVision FLEX Target Retrieval Solution, Low pH (50x) (Code K8005) for HIER.</p>
Product-specific limitations	Occasional labelling of tissue components in vessel walls and pancreatic stroma has been observed in immunohistochemistry.

Staining interpretation	Cells labeled by the antibody display a nuclear staining pattern except in mitotic cells, where the chromosomes and the cytoplasm are labeled.
Performance characteristics	<p><u>Normal tissues:</u> In small intestine and colon, mucosal neck cells and deep crypt cells are positive with the antibody. The same is the case with germinal centre cells of Payer's patches. Surface epithelium cells and other mucosal cells, e.g. Paneth's cells and cells of submucosal connective tissue are completely negative. The same is true for muscle cells, with the exception of a few positive cells in smooth muscles. Kidney, liver, pancreas and brain are negative (1).</p> <p><u>Abnormal tissues:</u> The antibody was in some cases found to stain membranes and cytoplasm. Membranous/cytoplasmic staining was reported in a study of 322 invasive ductal breast carcinomas (7). 24/24 prostate carcinomas were labeled by the antibody (8).</p>

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
References/ Bibliographie/ Literaturnachweise

1. Cattoretti G, Becker MHG, Key G, Duchrow M, Schlüter C, Galle J, et al. Monoclonal antibodies against recombinant parts of the Ki-67 antigen (MIB 1 and MIB 3) detect proliferating cells in microwave-processed formalin-fixed paraffin sections. J Pathol 1992;168:357-63.
2. Gerdes J, Lemke H, Baisch H, Wacker H-H, Schwab U, Stein H. Cell cycle analysis of a cell proliferation-associated human nuclear antigen defined by the monoclonal antibody Ki-67. J Immunol 1984;133:1710-5.
3. Gerdes J, Li L, Schlueter C, Duchrow M, Wohlenberg C, Gerlach C, et al. Immunobiochemical and molecular biologic characterization of the cell proliferation-associated nuclear antigen that is defined by monoclonal antibody Ki-67. Am J Pathol 1991;138:867-73.
4. Scholzen T, Gerdes J. The Ki-67 protein: from the known and the unknown [review]. J Cell Physiol 2000;182:311-22.
5. Key G, Kubbutat MH, Gerdes J. Assessment of cell proliferation by means of an enzyme-linked immunosorbent assay based on the detection of the Ki-67 protein. J Immunol Methods 1994;177:113-7.
6. Key G, Becker MHG, Baron B, Duchrow M, Schlüter C, Flad H-D, et al. New Ki-67-equivalent murine monoclonal antibodies (MIB 1-3) generated against bacterially expressed parts of the Ki-67 cDNA containing three 62 base pair repetitive elements encoding for the Ki-67 epitope. Lab Invest 1993;68:629-36.
7. Faratian D, Munro A, Twelves C, Bartlett JMS. Membranous and cytoplasmic staining of Ki67 is associated with HER2 and ER status in invasive breast carcinoma. Histopathology 2009;54:254-257.
8. Claudio PP, Zamparelli A, Garcia FU, Claudio L, Ammirati G, Farina A, et al. Expression of cell-cycle-regulated proteins pRb2/p130, p107, p27kip1, p53, mdm-2, and Ki-67 (MIB-1) in prostatic gland adenocarcinoma. Clin Cancer Res 2002; 8:1808-15.

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