**Monoclonal Mouse
Anti-Human
Collagen IV**
Clone CIV 22
Code M0785

Performance characteristics Normal tissues: The antibody shows the characteristic labeling of basement membranes in a variety of tissues and organs tested, including kidney, skin, striated and smooth muscle, spleen, lymph node, lung, placenta and tendon. In spleen and lymph nodes, the expected fragmented labeling of the discontinuous basement membranes of the sinusoids is observed, whereas other blood vessels exhibit a linear, continuous labeling. In kidneys, basement membranes of capillaries, parts of the mesangial matrix and the Bowman's capsule, and the tubular basement membranes are labeled by the antibody. The only basement membrane showing no labeling with the antibody is that of the corneal epithelium. All structures other than basement membranes are consistently not labeled with this antibody (2).

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Intended use

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human Collagen IV, Clone CIV 22, is intended for use in immunohistochemistry (IHC). The antibody labels type IV collagen and is a useful aid for the identification of basement membranes (1, 2). Differential classification of tumors is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

Summary and explanation

Collagen IV is a major constituent in particular of the lamina densa of basement membranes and ultrastructurally it looks amorphous. Basement membranes are thin extracellular matrices separating parenchymal, endothelial, and epithelial cells from underlying connective tissue. In glomerular and tubular basement membranes 40% of the protein is type IV collagen.

Detection of intracellular type IV collagen is generally difficult owing to a low concentration. In newly formed capillaries found in the inflammatory sites of rheumatoid arthritis synovium, intracellular type IV collagen may however be detected (1).

Structurally the protein consists of four domains. One of these, the triple helical collagen IV domain, is formed by the association of two polypeptide chains of type $\alpha 1$ and one of type $\alpha 2$. This domain is 340 nm long and is highly crosslinked by disulphide bridges (1, 2).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure, Materials Required, Not Supplied, Storage, Specimen Preparation, Staining Procedure, Quality Control, Troubleshooting, Interpretation of Staining, General Limitations.

Reagent provided

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 50 mmol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN₃.

Clone: CIV 22 (2). Isotype: IgG1, kappa.

Mouse IgG concentration: see label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

Purified pepsin fragments of type IV collagen isolated from human kidney (2).

Specificity

The specificity of the antibody has been assessed in radioimmunoassay (RIA) where it recognizes an epitope present in the human type IV collagen in native conformation, whereas it does not recognize the reduced and alkylated protein in the denatured state. Staining of immunoblots with the antibody was negative, further indicating that the antibody recognizes a conformational epitope on type IV collagen (2).

In immunoblotting or RIA, no cross-reactivity of the antibody with isolated human collagen types I, II, III and V could be detected (2).

Precautions

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin sections: The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labelling acetone fixed, frozen sections (3). The user must validate the staining procedure.

Staining procedure

These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.

Dilution: Monoclonal Mouse Anti-Human Collagen IV, Code M0785, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human kidney and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.

Visualization: Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.

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