**Polyclonal Rabbit
Anti-Human CD3**
Code A0452
ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Polyclonal Rabbit Anti-Human CD3, is intended for use in immunohistochemistry (IHC). The antibody labels CD3 and is a useful tool for the identification of T cells. Results aid in the classification of T-cell neoplasms (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Synonyms for antigen	T3, CD3 complex (3).
Summary and explanation	<p>CD3 consists of at least four different components ($\gamma, \delta, \epsilon, \zeta$) of 20-28 kDa. On the lymphocyte cell surface, CD3 is non-covalently associated with the T-cell receptor (TCR). It is believed that the CD3 components of the TCR/CD3 complex mediate signal transduction upon antigen recognition by TCR (4).</p> <p>CD3 is expressed by T cells in thymus, bone marrow, peripheral lymphoid tissue and blood (5, 6). The only other normal cell type known to be labelled by antibodies to CD3, is Purkinje cells in the cerebellum (7). The CD3 antigen is first detectable in early thymocytes and its appearance probably represents one of the earliest signs of commitment to the T-cell lineage (5). In immature thymocytes the CD3 expression is solely cytoplasmic, membrane-associated CD3 appears subsequently upon cell maturation (5).</p> <p>The majority of T-cell neoplasms may express the CD3 antigen, while it appears to be absent from non-T-cell lymphoid malignancies (8). Consistent with the pattern of synthesis of the antigen in normal thymocytes, the earliest site where CD3 is detectable within neoplastic cells is in the cell cytoplasm (5).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Affinity-isolated rabbit anti-human CD3 provided in liquid form. In 0.05 mol/L Tris-HCl, 0.1 mol/L NaCl, 15 mmol/L NaN₃. The affinity isolation has been performed using immobilized CD3 peptide.</p> <p><u>Protein concentration:</u> See label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Synthetic peptide comprising amino acids 156-168 from the cytoplasmic part of the human CD3 ϵ -chain coupled to thyroglobulin (1).
Specificity	<p>In Western blotting, the antibody detects bands of the expected molecular weights for CD3 antigens (2). The antibody recognizes CD3ϵ in both a T-cell line (Jurkat) and a natural killer cell line (NK11), but does not react with lysates prepared from several B-cell lines (Raji, Ramos and JY), a myeloid cell line (U937) or a colon carcinoma cell line (Colo-205) (9).</p> <p>In immunoprecipitation from Nonidet P40 lysates of surface-iodated T lymphoblasts, the antibody precipitates the γ (26 kDa), δ (21 kDa) and ϵ (19 kDa) chain of the CD3 molecule, similar to the precipitation pattern seen using the well-characterized monoclonal mouse anti-human CD3, clone UCHT1 (1).</p> <p>In ELISA, the antibody labels the CD3 peptide used as immunogen.</p>
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, pH 9, Code S2368. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or Dako Target Retrieval Solution, Citrate pH 6, Code S2369, and pre-treatment with Dako Proteinase K, Code S3020. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling frozen sections (1). The user must validate the staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Polyclonal Rabbit Anti-Human CD3, Code A0452, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, pH 9, Code S2368, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to the same protein concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen.</p>

Visualization: Dako EnVision+HRP kits, e.g. Code K4009, are recommended. Follow the procedure enclosed with the selected visualization kit.

Staining interpretation

Cells labeled by the antibody display staining confined to the cytoplasm and/or cell membrane.

Performance characteristics

Normal tissues: The antibody labels T cells in various tissues including tonsil and colon.

Abnormal tissues: The antibody labeled 73/96 T-cell neoplasms, including 7/9 lymphoblastic, 25/35 pleomorphic, 5/5 immunoblastic, 5/5 angioimmunoblastic lymphadenopathy-like, 2/2 T-zone lymphomas, 19/19 mycosis fungoides/Sézary syndrome, 2/3 Lennert's lymphomas, 4/13 Ki-1 positive anaplastic large cell lymphomas, 3/4 lymphomatoid papulosis, 1/1 coeliac-associated lymphoma (1). The antibody labeled 149/149 cases of T-cell (precursor) acute lymphoblastic leukaemias/lymphomas (ALL). In 131/149 cases, 100% of the cells were labeled; in 14 cases between 50-90% of the cells were labeled; and in 4 cases less than 50% of the cells were labeled. No labeling was observed in 68/68 cases of B-cell (precursor) ALL, apart from reactive infiltrating T cells (2).

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

- Mason DY, Cordell J, Brown M, Pallesen G, Ralfkiaer E, Rothbard J, et al. Detection of cells in paraffin wax embedded tissue using antibodies against a peptide sequence from the CD3 antigen. J Clin Pathol 1989;42:1194-1200.
- Pilozzi E, Pulford K, Jones M, Müller-Hermelink H-K, Falini B, Ralfkiaer E, et al. Co-expression of CD79a (JCB117) and CD3 by lymphoblastic lymphoma. J Pathol 1988;186:140-3.
- CD guide In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1111.
- Jacobs H. Pre-TCR/CD3 and TCR/CD3 complexes: decamers with differential signalling properties? Immunol Today 1997;18:565-9.
- Campana D, Thompson JS, Amlot P, Brown S, Janossy G. The cytoplasmic expression of CD3 antigens in normal and malignant cells of the T lymphoid lineage. J Immunol 1987;138:648-55.
- Tunnacliffe HA, Olsson C, Traunecker A, Krissansen GW, Karjalainen K, De la Hera A. The majority of CD3 epitopes are conferred by the epsilon chain. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 295-6.
- Garson JA, Beverley PC, Coakham HB, Harper EI. Monoclonal antibodies against human T lymphocytes label Purkinje neurones of many species. Nature 1982;298:375-7.
- Erber WN, Mynheer LC, Mason DY. APAAP labeling of blood and bone-marrow samples for phenotyping leukaemia. Lancet 1986;i:761-5.
- Lanier LL, Chang C, Spitz H, Phillips JH. Expression of cytoplasmic CD3ε proteins in activated human adult natural killer (NK) cells and CD3γ, δ, ε complexes in fetal NK cells. Implications for the relationship of NK and T lymphocytes. J Immunol 1992;149:1876-80.

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