Code M0755

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD20cy, Clone L26, is intended for use in immunohistochemistry (IHC). The antibody labels cells of the B-cell lineage and is a useful aid for the classification of neoplasms of B-cell derivation (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Synonyms for antigen	L26 (2).
Summary and explanation	<p>CD20 is a transmembrane, non-glycosylated protein expressed on B-cell precursors and mature B cells, but is lost following differentiation into plasma cells (3). In resting B cells, CD20 appears in a 33 kDa non-phosphorylated form. After mitogen stimulation, CD20 becomes heavily phosphorylated (35-37 kDa isoforms), and it is a dominant phosphoprotein in activated B cells, B-cell lines, and hairy cell leukemias (2). The long N- and C-terminal ends of the protein are located on the cytoplasmic side of the membrane and only a minor portion of the protein is exposed on the cell surface (3). Antibodies reacting with CD20 cytoplasmic epitopes are designated CD20cy (2). It is suggested that CD20 plays a direct role in regulating the transmembrane conductive Ca²⁺ flux of B cells which indicates a possible function for CD20 as a regulator of proliferation and differentiation (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure, Materials Required, Not Supplied, Storage, Specimen Preparation, Staining Procedure, Quality Control, Troubleshooting, Interpretation of Staining, General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L Na₃N.</p> <p><u>Clone:</u> L26 (1, 4). <u>Isotype:</u> IgG2a, kappa.</p> <p><u>Mouse Ig concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Human tonsil B cells (4).
Specificity	<p>The antibody was clustered as anti-CD20 at the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Boston 1993 (2).</p> <p>SDS-PAGE analysis of immunoprecipitates formed between ¹²⁵I-labeled tonsil cell lysate and the antibody shows reaction primarily with 30 kDa and 33 kDa polypeptides (4).</p> <p>Studies using COS-1 cells transfected with cDNA encoding the CD20 molecule, indicate that the antibody labels an intracytoplasmic epitope localized on the CD20 molecule (5).</p>
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (Na₃N), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin or B5-fixative (6). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is recommended. For tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed frozen sections (1) and acetone-fixed cell preparations (1, 4). The user must validate the staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human CD20cy, Code M0755, may be used at a dilution range of 1:200-1:400 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG2a, Code X0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Visualization:</u> Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p>
Product-specific limitations	Rare cases of CD20-expressing peripheral T-cell lymphomas have been reported (7, 8).

Staining interpretation	B cells labeled by the antibody display staining of the cytoplasmic side of the cell surface membrane.
Performance characteristics	<p>Normal tissues: In normal lymphoid tissue, the antibody labeled germinal centre cells, mantle zone lymphocytes, and scattered interfollicular lymphocytes, but not T cells, histiocytes and plasma cells (6, 8). No labeling was observed in epidermis, sebaceous glands, hair follicles and eccrine glands in the skin, follicular epithelium in the thyroid, pneumocytes and bronchial epithelium of the lung, and a large number of other normal non-lymphoid tissues tested (6).</p> <p>Abnormal tissues: Reaction with the antibody was revealed in most of 131 B-cell neoplasms tested (1). Labeling with the antibody shows that in the differentiation of B cells, the CD20 antigen is not expressed on very immature lymphoid cells (0/6 acute undifferentiated leukemias), but begins to be expressed on early maturational stages (14/34 common acute lymphoblastic and 7/9 pre-B acute lymphoblastic leukemias), and then, the CD20 antigen is fully expressed on mature B cells (15/15 chronic lymphocytic, 3/3 prolymphocytic, 3/3 hairy cell, 6/7 lymphosarcoma cell leukemias, and 45/46 B-cell malignant lymphomas including Burkitt, Waldenström, and immunoblastic B-cell lymphomas. The CD20 antigen disappears on plasma cells, and only 1/2 plasma cell leukemias and 0/12 myelomas were labeled by the antibody (1). Other studies have provided comparable results showing labeling of 44/44 large cell and immunoblastic B-cell lymphomas (6), and all of 40 B-cell lymphomas, with the exception of common acute lymphoblastic leukemias and malignant lymphoma plasmacytic. In Hodgkin's disease, strong surface membrane staining of Reed-Sternberg cells was observed in 9/27 cases (7). Of lymphoproliferative diseases of T-cell lineage, 0/73 were labeled by the antibody (1), whereas other studies showed 1/18 (7) and 1/111 (8) labeled cases.</p>

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
References / Bibliographie / Literaturnachweise

1. Takami T, Qi C-F, Yamada T, Yamashina M, Kon S-I, Ishii Y, et al. B20.3. Reactivity and specificity of L26 (pan-B-cell mAb) on 322 cases of fresh and paraffin-embedded lymphoproliferative diseases. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 134-6.
2. Zhou L-J, Tedder TF. CD20 workshop panel report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 511-4.
3. Tedder TF, Engel P. CD20: a regulator of cell-cycle progression of B lymphocytes. Immunology Today 1994;15:450-4.
4. Ishii Y, Takami T, Yuasa H, Takei T, Kikuchi K. Two distinct antigen systems in human B lymphocytes: identification of cell surface and intracellular antigens using monoclonal antibodies. Clin Exp Immunol 1984;58:183-92.
5. Mason DY, Comans-Bitter WM, Cordell JL, Verhoeven MAJ, van Dongen JJM. Antibody L26 recognizes an intracellular epitope on the B-cell-associated CD20 antigen. Am J Pathol 1990;136:1215-22.
6. Cartun RW, Coles FB, Pastuszak WT. Utilization of monoclonal antibody L26 in the identification and confirmation of B-cell lymphomas. A sensitive and specific marker applicable to formalin- and B5-fixed, paraffin-embedded tissues. Am J Pathol 1987;129:415-21.
7. Norton AJ, Isaacson PG. Monoclonal antibody L26: an antibody that is reactive with normal and neoplastic B lymphocytes in routinely fixed and paraffin wax embedded tissue. J Clin Pathol 1987;40:1405-12.
8. Blakolmer K, Vesely M, Kummer JA, Jurecka W, Mannhalter C, Chott A. Immunoreactivity of B-cell markers (CD79a, L26) in rare cases of extranodal cytotoxic peripheral T- (NK/T-) cell lymphomas. Mod Pathol 2000;13:766-72.

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