**Monoclonal Mouse  
Anti-Human  
CD44, Phagocytic Glycoprotein-1**  
Clone DF1485  
Code M7082

**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use.  Monoclonal Mouse Anti-Human CD44, Phagocytic Glycoprotein-1, Clone DF1485, is intended for use in immunohistochemistry (IHC). The antibody labels lymphocytes, both T and B cells, monocytes, granulocytes, erythrocytes, epithelial cells (1) and mast cells (2). Differential classification of tumors is aided by a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains
<b>Synonyms for antigen</b>	CD44s, ECMRIII, Hermes antigen, homing-associated cell adhesion molecule (H-CAM), Pgp-1 (2-4).
<b>Summary and explanation</b>	CD44 is a type I transmembrane adhesion molecule with a molecular mass of 85-90 kDa. It is a glycoprotein with extensive O-linked glycosylation and 6 potential N-linked glycosylation sites. A portion of the N-terminal domain of CD44 shows homology with a part of cartilage proteoglycan core and link proteins (5). CD44 is a cell surface receptor for hyaluronate suggesting a role in the regulation of cell and cell-substrate interactions as well as cell migration (3).  CD44 is expressed on a wide variety of cells. In the hematopoietic system it is present on monocytes, granulocytes, erythrocytes, B cells and mature T cells (1, 2). Outside the hematopoietic system it is expressed on epithelial cells, fibroblasts, skeletal muscle, central nervous system white matter, and a wide variety of tumors (5, 6).  Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure, Materials Required, Not Supplied, Storage, Specimen Preparation, Staining Procedure, Quality Control, Troubleshooting, Interpretation of Staining, General Limitations.
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> , <u>Clone:</u> DF1485. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.  The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
<b>Specificity</b>	Anti-CD44, DF1485, was included in the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens (1), and a number of laboratories have confirmed the reactivity against CD44.
<b>Precautions</b>	1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700. Less optimal results are obtained with Dako Target Retrieval Solution, citrate pH 6, Code S2369 or Dako Target Retrieval Solution, pH 9, Code S2368, while pre-treatment of tissues with Dako Proteinase K, Code S3020, was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.
<b>Staining procedure</b>	These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.  <u>Dilution:</u> Monoclonal Mouse Anti-Human CD44, Phagocytic Glycoprotein-1, Code M7082, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG1 concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S 0809.  <u>Quality control:</u> Positive and negative tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen.  <u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

<b>Staining interpretation</b>	Cells labeled by the antibody display staining confined to the cell surface membrane (6). However, cytoplasmic staining has been reported (4, 7).
<b>Performance characteristics</b>	<p><b>Normal tissues:</b> The antibody labels approximately 90% of all lymphocytes, both T cells and B cells. It also labels monocytes, granulocytes, erythrocytes and epithelial cells. The epitope is weakly expressed on platelets (1, 6).</p> <p><b>Abnormal tissues:</b> The antibody has been shown to label both normal and neoplastic tissue mast cells (2). The expression of CD44 was lower in carcinomas than adenomas of the colorectum and appendix, respectively. The antibody labeled 68/71 (97%) of colorectal adenomas and 34/39 (87%) of appendiceal adenomas, whereas the antibody labeled 109/137 (80%) of colorectal carcinomas and 19/34 (56%) of appendiceal carcinomas (4). In a study of 49 cases of colorectal carcinoma labeling was found in 41% of malignant epithelial tissue as compared to 10% of normal tissue. The expression of CD44 in the epithelium was higher in rectal cancer (48%) than in colon cancer (21%) (7).</p>

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

1. Stoll M, Dalchau R, Schmidt RE. N6 Cluster report: CD44. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 619-22.
2. Horny H-P, Menke DM, Kaiserling E. Neoplastic human tissue mast cells express the adhesion molecule CD44/HCAM. Virchows Arch 1996;429:91-4.
3. Aruffo A, Stamenkovic I, Melnick M, Underhill CB, Seed B. CD44 is the principal cell surface receptor for hyaluronate. Cell 1990;61:1303-13.
4. Carr NJ, Emory TS, Sobin LH. Epithelial neoplasms of the appendix and colorectum. An analysis of cell proliferation, apoptosis, and expression of p53, CD44, and bcl-2. Arch Pathol Lab Med. 2002;126:837-41.
5. Stamenkovic I, Amiot M, Pesando JM, Seed B. A lymphocyte molecule implicated in lymph node homing is a member of the cartilage link protein family. Cell 1989; 56:1057-62.
6. Dalchau R, Flanagan BF, Allen AK, Daar AS, Fabre JW. Chemical composition and tissue distribution of the human CD44 glycoprotein. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 622-5.
7. Khoursheed M, Mathew TC, Makar RR, Sonia L, Abul H, Asfar S, et al. Expression of CD44s in human colorectal cancer. Pathol Oncol Res. 2002;8:170-4.

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