**Negative Control
Rabbit Immunoglobulin Fraction
(Solid-Phase Absorbed)
Code X0936**

ENGLISH	
Intended use	For in vitro diagnostic use. Negative Control, Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed) is intended for use as negative control for Dako solid-phase absorbed polyclonal primary antibodies in immunohistochemistry, immunoblotting and other techniques. The reagent has been produced according to the same method as Dako solid-phase absorbed polyclonal primary antibodies, thus, it is well-suited as a negative control for these antibodies.
Reagent provided	Purified immunoglobulin fraction of serum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN ₃ . The immunoglobulin (Ig) fraction is isolated from serum of healthy, non-immunized rabbits (Danish Whites) by salting out and ion exchange chromatography. The Ig fraction is solid-phase absorbed with human plasma. This reduces the non-specific background, and makes the reagent a useful control for solid-phase absorbed antibodies. <u>Protein concentration g/L:</u> See label on vial.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. The product may be used in different techniques and in combination with different sample types and materials, therefore each individual laboratory should validate the test system applied. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Dilution guidelines	X0936 should be diluted to the same protein concentration as the test antibody. Unless the stability of the diluted negative control has been established in the actual staining procedure, it is recommended to dilute the reagent immediately before use.

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