**Monoclonal Mouse  
Anti-Human CD8**  
Clone C8/144B

**Code M7103**

**ENGLISH**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD8, Clone C8/144B, is intended for use in immunohistochemistry (IHC). The antibody labels cytotoxic/suppressor T cells. Results aid in the classification T-cell lymphoma (1) and mycosis fungoides (2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
<b>Summary and explanation</b>	<p>CD8 is a transmembrane glycoprotein with a molecular mass of 68 kDa. It is expressed as a disulfide-linked heterodimer comprising a 32-34 kDa <math>\alpha</math>- and a 30-32 kDa <math>\beta</math>-chain, or as a homodimer comprising two <math>\alpha</math>-chains. Both CD8<math>\alpha</math> and CD8<math>\beta</math> have a typical immunoglobulin variable region-like domain in an N-terminal extracellular portion that makes them members of the immunoglobulin gene superfamily. CD8 is expressed mostly as the <math>\alpha\beta</math> heterodimer by a majority of thymocytes, and by class I major histocompatibility complex restricted, mature, suppressor/cytotoxic T cells. A proportion of <math>\gamma\delta</math> T cells and NK cells express the CD8 <math>\alpha\alpha</math> homodimer (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
<b>Reagent provided</b>	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN<sub>3</sub>.</p> <p><u>Clone:</u> C8/144B (1). <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
<b>Immunogen</b>	<p>Synthetic peptide corresponding to the 13 C-terminal amino acids of cytoplasmic domain of human CD8<math>\alpha</math> coupled to thyroglobulin (1).</p>
<b>Specificity</b>	<p>SDS-PAGE analysis of immunoprecipitates formed between lysates of <sup>125</sup>I-labeled human T lymphoblasts and the antibody shows reaction primarily with a 32 kDa polypeptide corresponding to CD8<math>\alpha</math> (1).</p>
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For in vitro diagnostic use.</li> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> <li>Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
<b>Specimen preparation</b>	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin or Bouin's fixative (1). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed, frozen sections or cell preparations (1, 2). The user must validate the staining procedure.</p>
<b>Staining procedure</b>	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human CD8, Code M7103, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen.</p> <p><u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>

<b>Product-specific limitations</b>	The antibody labels native cytokeratin 15 in human hair follicle bulge keratinocytes and a discrete area of the outer root sheet basal layer near the attachment of the arrector pili muscle without labeling the remaining hair follicle (4).
<b>Staining interpretation</b>	Cells labeled by the antibody display staining confined to the surface membrane.
<b>Performance characteristics</b>	<u>Normal tissues:</u> The antibody labels cytotoxic/suppressor T cells in T-cell areas of spleen and tonsil, and splenic sinusoidal lining cells (1). In colorectal biopsies, the antibody labels cytotoxic/suppressor T cells in rectum, sigmoid-, descending-, transverse-, and ascending colon, and caecum (5). <u>Abnormal tissues:</u> In mycosis fungoides, the CD8 antibody together with anti-CD4 demonstrated an elevated CD4:CD8 ratio in the epidermis compared with the ratio in inflammatory lesions (2).

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
1. Mason DY, Cordell JL, Gaulard P, Tse AGD, Brown MH. Immunohistological detection of human cytotoxic/suppressor T cells using antibodies to a CD8 peptide sequence. J Clin Pathol 1992;45:1084-8.
2. Nuckols JD, Shea CR, Horenstein MG, Burchette JL, Prieto VG. Quantitation of intraepidermal T-cell subsets in formalin-fixed, paraffin-embedded tissue helps in the diagnosis of mycosis fungoides. J Cutan Pathol 1999;26:169-75.
3. Nakauchi H. TC9. CD8 Workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 65-7.
4. Lyle S, Christofidou-Solomidou M, Liu Y, Elder DE, Albelda S, Cotsarelis G. Human hair follicle bulge cells are biochemically distinct and possess an epithelial stem cell phenotype. J Invest Dermatol Symp Proc 1999;4:296-301.
5. Yamagata K, Tanaka M, Kudo H. A quantitative immunohistochemical evaluation of inflammatory cells at the affected and unaffected sites of inflammatory bowel disease. J Gastroenterol 1998;13:801-8.

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