



## Certificate of Analysis

### Anti-p75NTR (Neurotrophin Receptor), clone ME20.4

(mouse monoclonal IgG<sub>1</sub>)

Catalog # 05-446

Lot # 2484171

**Immunogen:** Human melanoma cell line WM245.  
Hybridoma Clone ME20.4

**Specificity:** Recognizes the low affinity p75NTR,  
Mr 69kDa.

**Species Cross-reactivity:** Human and monkey.

**Formulation:** 100µg of protein G purified mouse IgG<sub>1</sub>  
in 132µl of 0.1M Tris-glycine, pH 7.4, 0.15M NaCl,  
0.05% sodium azide before the addition of glycerol to  
30%. Liquid at -20°C.

**Storage and Stability:** Stable for 2 years at -20°C  
from date of shipment. For maximum recovery of  
product, centrifuge the original vial after thawing and  
prior to removing the cap.

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC  
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

#### Quality Control Testing


**Immunoblot Analysis:** 0.1-1µg/ml of this lot  
detected p75NTR in RIPA lysates from HS294T  
cells.

**Immunoprecipitation:** 4µg of a previous lot immuno-  
precipitated p75NTR from 500µg of HS294T RIPA  
lysate.

#### Additional Research Applications

**Immunohistochemistry:** Reported to detect p75NTR in  
tissue section. Also used in immuno-electron  
microscopy.

**Biological Blockade:** Reported to block <sup>125</sup>I-labeled NGF  
binding to its receptor.



#### Immunoblot Analysis

Representative blot from a previous lot.  
HS294T cell lysate was resolved by  
electrophoresis, transferred to nitrocellulose  
and probed with anti-p75NTR (1µg/ml).  
Proteins were visualized using a goat anti-  
mouse secondary antibody conjugated to  
HRP and a chemiluminescence detection  
system. Arrow indicates p75 NTR (69kDa).

#### Application References:

1. Ross, A.H., *et al.*, *Proc. Natl. Acad. Sci. USA* **81**: 6681-6685, 1984.
2. Dominici, C., *J. Neurooncol.* **31**: 57-64, 1997.
3. Wakabayashi, Y., *Neurosci. Lett.* **186**: 9-12, 1995.
4. Caneva, L., *Blood Cells Mol. Dis.* **21**: 73-85, 1995.

### Immunoblot Protocol

1. Perform SDS-polyacrylamide gel electrophoresis (SDS-PAGE) on a cell lysate sample (cell lysis buffer: 50mM Tris-HCl, pH7.4; 1% NP-40; 0.25% sodium deoxycholate; 150mM NaCl; 1mM EDTA; 1mM PMSF; 1 $\mu$ g/ml each aprotinin, leupeptin, pepstatin; 1mM Na<sub>3</sub>VO<sub>4</sub>; 1mM NaF) and transfer the proteins to nitrocellulose. Wash the blotted nitrocellulose twice with water.
2. Block the blotted nitrocellulose in freshly prepared PBS containing 3% nonfat dry milk (Catalog # 20-200), (PBS-MLK) for 20 minutes at room temperature with constant agitation.
3. Incubate the nitrocellulose with **0.1-1 $\mu$ g/ml of anti-p75NTR (Neurotrophin Receptor)**, diluted in freshly prepared PBS-MLK overnight with agitation at 4°C.
4. Wash the nitrocellulose twice with water.
5. Incubate the nitrocellulose in the secondary reagent of choice (a **goat anti-mouse** HRP conjugated IgG, Catalog # 12-349, 1:2000 dilution was used) in PBS-MLK for 1.5 hours at room temperature with agitation.
6. Wash the nitrocellulose with water twice.
7. Wash the nitrocellulose in PBS-0.05% Tween 20 for 3-5 minutes.
8. Rinse the nitrocellulose in 4-5 changes of water.
9. Use detection method of choice (enhanced chemiluminescence was used).

### Immunoprecipitation Protocol

1. Dilute the cell lysate before beginning the immunoprecipitation to roughly 1 $\mu$ g/ $\mu$ l total cell protein in a microcentrifuge tube with PBS.
2. Add **4 $\mu$ g of anti-p75NTR (Neurotrophin Receptor)** to 500 $\mu$ g-1mg cell lysate.
3. Gently rock the reaction mixture at 4°C overnight.
4. Capture the immunocomplex by adding 100 $\mu$ l (50 $\mu$ l packed beads) of washed Protein G agarose bead slurry (Catalog # 16-266).
5. Gently rock the reaction mixture at 4°C for 2 hours.
6. Collect the agarose beads by pulsing (5 seconds in the microcentrifuge at 14,000 x g), and drain off the supernatant. Wash the beads 3 times with either ice-cold cell lysis buffer or PBS.
7. Resuspend the agarose beads in 60 $\mu$ l 2X Laemmli sample buffer.
8. The agarose beads can either be frozen for later use or boiled for 5 minutes. Collect the beads by a microcentrifuge pulse. SDS-PAGE and subsequent immunoblot analysis can be performed on a sample of the supernatant.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©1997 – 2013 Merck KGaA, Darmstadt. All rights reserved. No part of these works may be reproduced in any form without permission in writing.